

Headwalls

NuLOOK
Refresh




Amico

Refresh Your Room

IN UNDER
4 HOURS

Without **Breaking the Wall**
or **Cutting Pipes**



 Call Us for a Demonstration 1.877.462.6426

After



**See
How**

The NuLook Refresh Process



Original Hospital Room

We test the performance of your current system.



Remove All Existing Fixtures

Without breaking into the wall or cutting pipes.



Hang the Prefabricated Frame

The NuLook frame and adapters are then installed, retrofitting, and relocating your existing outlets.



Attach Panels and Complete Your NuLook Room

The entire process took less than 4 hours!



NuLook for Any Condition

INSTALL OVER EXISTING CONSOLES

Increase your room's functionality by installing over top of existing gas console boxes while recessing the bed location to maximize bed clearance. Integrate recessed equipment rails for better equipment management.



RETROFIT UNSERVICEABLE HEADWALLS

Replace headwalls made by past manufacturers who no longer provide parts or service. NuLook allows you to bring unserviceable headwalls back into operation in live areas, all without construction or recertifying.



NuLook for Every Need

BE UNIQUE, CUSTOMIZE YOUR DESIGN

Transform your LDR with our custom design capabilities, regardless of the existing conditions. Custom graphics, lighting, and layout allows you to make your rooms truly unique without exceeding your budget.



UPGRADE YOUR ROOM AND DESIGN


Increase electrical circuitry and create a home-like environment. Maximize patient and nurse comfort with graphic panels and integrated bedside cabinets.



NuLook for Every Acuity

INCREASE PATIENT AND NURSE SATISFACTION

With NuLook anything is possible; conceal gases behind sliding panels. Installs in any existing condition, increasing your patient and nurse satisfaction, while avoiding the time and cost of a traditional gut and renovation.



CONVERT EXISTING SEMI-PRIVATE ROOMS INTO PRIVATE

The system has been designed to install in any condition. We floor mount our unit when working with pre-existing plaster walls in order to provide greater structural support. NuLook gives you the flexibility to convert your rooms from semi-private into private.



NuLook for Every Budget

UPGRADE WITH LESS DOWNTIME

Improve infection control, usability, and aesthetics in your existing OB suites. NuLook's fully integral design process allows us to achieve your goals without the need to recertify.




LIMITLESS DESIGN CAPABILITIES

NuLook comes in all shapes and sizes. Update your outlets, increase infection control, and stay within your budget.



The Possibilities are Endless





Bringing Non-Disruptive Change to Your Facility

Refresh

Construction-free headwall system that installs without the need to recertify, allowing for non-disruptive change to your existing facility.

Casework

Future adaptable and demountable casework prefabricated with industry leading durability, warranty, and flexibility.

FlexWall

Acuity-adaptable prefabricated headwall system that allows for future flexibility, while simplifying the healthcare construction process.

Skeleton

Prefabricated stud-frame system engineered to streamline the conventional headwall construction process.

Accessories

Accessories

Price Book 2020



Basket 5 x 5 x 4

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

5" (12.7cm) x 5" (12.7cm) x 4" (10.2cm) (w x d x h)

Small storage basket great for BP cuffs and miscellaneous odds and ends. Designed to sit below the horizontal rail, or in the center for vertical rail systems.



050504-(X)-BB

A - Amico Horizontal Rail	D - Ohmeda Dovetail Vertical Rail	E25 - DIN 25mm Rail	R200 - 2" Diam. Pole Clamp
F - Fairfield Rail	G - Hill-Rom Vertical Rail	DW - Direct To Wall	MON - Vertical Monitor Channel
V - Amico Vertical Rail	S - HSI Vertical Rail	R125 - 1.25" Diam. Pole Clamp	

Part#

050504-A-BB
050504-F-BB
050504-V-BB
050504-D-BB
050504-G-BB
* 050504-S-BB

Part#

050504-E25-BB
050504-DW
§ 050504-R125-BB
050504-R200-BB
050504-MON-BB

§ Comes with cord wrap on the side of the basket.

* Long lead time product.

Fasteners required to mount to wall are not provided.

Basket 9 x 9 x 10

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

9" (22.9cm) x 9" (22.9cm) x 10" (25.4cm) (w x d x h)

Medium storage basket designed to accommodate small objects for nurses and doctors convenience. Designed to sit below the horizontal rail, or in the center for vertical rail systems.



090910-(X)-BB

A - Amico Horizontal Rail	D - Ohmeda Dovetail Vertical Rail	E25 - DIN 25mm Rail	R200 - 2" Diam. Pole Clamp
F - Fairfield Rail	G - Hill-Rom Vertical Rail	DW - Direct To Wall	MON - Vertical Monitor Channel
V - Amico Vertical Rail	S - HSI Vertical Rail	R125 - 1.25" Diam. Pole Clamp	

Part#

090910-A-BB
090910-F-BB
090910-V-BB
090910-D-BB
090910-G-BB
* 090910-S-BB

Part#

090910-E25-BB
090910-DW
§ 090910-R125-BB
090910-R200-BB
090910-MON-BB

§ Comes with cord wrap on the side of the basket.

* Long lead time product.

Fasteners required to mount to wall are not provided.

Basket 10 x 6 x 9

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

10" ^(25.4cm) x 6" ^(15.2cm) x 9" ^(22.9cm) (w x d x h)

Medium basket that works well with small items in patient rooms can also work as a convenient waste receptacle. Designed to sit below the horizontal rail, or in the center for vertical rail systems.



100609-(X)-BB

A - Amico Horizontal Rail	D - Ohmeda Dovetail Vertical Rail	E25 - DIN 25mm Rail	R200 - 2" Diam. Pole Clamp
F - Fairfield Rail	G - Hill-Rom Vertical Rail	DW - Direct To Wall	MON - Vertical Monitor Channel
V - Amico Vertical Rail	S - HSI Vertical Rail	R125 - 1.25" Diam. Pole Clamp	

Part#

100609-A-BB
100609-F-BB
100609-V-BB
100609-D-BB
100609-G-BB
* 100609-S-BB

Part#

100609-E25-BB
100609-DW
§ 100609-R125-BB
§ 100609-R200-BB
100609-MON-BB

§ Comes with cord wrap on the side of the basket.

* Long lead time product.

Fasteners required to mount to wall are not provided.

Basket 10 x 10 x 4

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

10" ^(25.4cm) x 10" ^(25.4cm) x 4" ^(10.2cm) (w x d x h)

Great for miscellaneous storage and organizing supplies. Designed to sit below the horizontal rail, or in the center for vertical rail systems.



101004-(X)-BB

A - Amico Horizontal Rail	D - Ohmeda Dovetail Vertical Rail	E25 - DIN 25mm Rail	R200 - 2" Diam. Pole Clamp
F - Fairfield Rail	G - Hill-Rom Vertical Rail	DW - Direct To Wall	MON - Vertical Monitor Channel
V - Amico Vertical Rail	S - HSI Vertical Rail	R125 - 1.25" Diam. Pole Clamp	

Part#

101004-A-BB
101004-F-BB
101004-V-BB
101004-D-BB
101004-G-BB
* 101004-S-BB

Part#

101004-E25-BB
101004-DW
§ 101004-R125-BB
101004-R200-BB
101004-MON-BB

§ Comes with cord wrap on the side of the basket.

* Long lead time product.

Fasteners required to mount to wall are not provided.

Basket 10 x 10 x 7

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

101007-(X)-BB

10" _(25.4cm) x 10" _(25.4cm) x 7" _(17.8cm) (w x d x h)

Medium basket to keep medical products organized and within reach. Designed to sit below the horizontal rail, or in the center for vertical rail systems.



A - Amico Horizontal Rail	D - Ohmeda Dovetail Vertical Rail	E25 - DIN 25mm Rail	R200 - 2" Diam. Pole Clamp
F - Fairfield Rail	G - Hill-Rom Vertical Rail	DW - Direct To Wall	MON - Vertical Monitor Channel
V - Amico Vertical Rail	S - HSI Vertical Rail	R125 - 1.25" Diam. Pole Clamp	

Part#

101007-A-BB

101007-F-BB

101007-V-BB

101007-D-BB

101007-G-BB

* 101007-S-BB

Part#

101007-E25-BB

101007-DW

§ 101007-R125-BB

101007-R200-BB

101007-MON-BB

§ Comes with cord wrap on the side of the basket.

* Long lead time product.

Fasteners required to mount to wall are not provided.

Basket 13 x 9 x 7

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

130907-(X)-BB

13" _(33.0cm) x 9" _(22.9cm) x 7" _(17.8cm) (w x d x h)

Easily accessible, and versatile storage basket which can be effortlessly transferred around the room or removed for sterilization. Designed to sit below the horizontal rail, or span across two vertical rails.



A - Amico Horizontal Rail	E25 - DIN 25mm Rail	MON - Vertical Monitor Channel	V138 - 14" Amico Recessed Flatwall Spanning 2 Vertical Rails
F - Fairfield Rail	DW - Direct To Wall	R200 - 2" Diam. Pole Clamp	

Part#

130907-A-BB

130907-F-BB

130907-E25-BB

130907-R200-BB

Part#

130907-DW

130907-MON-BB

^ 130907-V138-BB

Fasteners required to mount to wall are not provided.

^ 13.88" (35.2cm) between centers of two mounting points. Not compatible with Modular.

Basket 18 x 9 x 7

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

18" (45.7cm) x 9" (22.9cm) x 7" (17.8cm) (w x d x h)

Large storage basket useful for those hard to store items that take up space on a shelf or in a cupboard. Designed to sit below the horizontal rail.



180907-(X)-BB

A - Amico Horizontal Rail

F - Fairfield Rail

E25 - DIN 25mm Rail

DW - Direct To Wall

Part#

180907-A-BB
180907-F-BB

Part#

180907-E25-BB
180907-DW

Fasteners required to mount to wall are not provided.

Basket 19 x 5 x 3

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

19" (48.3cm) x 5" (12.7cm) x 3" (7.6cm) (w x d x h)

Storage basket that accommodates larger items to free up space throughout the room. Designed to sit below the horizontal rail, or span across two vertical rails.



190503-(X)-BB

A - Amico Horizontal Rail

F - Fairfield Rail

E25 - DIN 25mm Rail

DW - Direct To Wall

V138 - 14" Amico Recessed Flatwall
Spanning 2 Vertical Rails

V178 - 19" Amico Surface Flatwall
Spanning 2 Vertical Rails

V188 - 19" Amico Recessed Flatwall
Spanning 2 Vertical Rails

Part#

190503-A-BB
190503-F-BB
190503-E25-BB

190503-DW

Part#

^ 190503-V138-BB
^ 190503-V178-BB
^ 190503-V188-BB

Fasteners required to mount to wall are not provided.

^ 13.88"(35.2cm) between centers of two mounting points for 190503-V138-BB. Not compatible with Modular.

17.81"(45.2cm) between centers of two mounting points for 190503-V178-BB. Not compatible with Modular.

18.88"(47.9cm) between centers of two mounting points for 190503-V188-BB. Not compatible with Modular.

Basket 24 x 5 x 3

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

24"^(61.0cm) x 5"^(12.7cm) x 3"^(7.6cm) (w x d x h)

Large storage basket that can be moved and relocated easily on any part of the rail. Designed to sit below the horizontal rail.



240503-(X)-BB

A - Amico Horizontal Rail
F - Fairfield Rail

E25 - DIN 25mm Rail

DW - Direct To Wall
E25 - DIN 25mm Rail

V223 - 24" Surface Flatwall Spanning 2
Vertical Rails

Part#

240503-A-BB
240503-F-BB
240503-E25-BB

Part#

240503-DW
^ 240503-V223-BB

^ 22.3" (56.6cm) between centers of two mounting points. Not compatible with Modular.

Fasteners required to mount to wall are not provided.

Basket 24 x 9 x 7

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

24"^(61.0cm) x 9"^(22.9cm) x 7"^(17.8cm) (w x d x h)

Storage basket with ample room for larger items that require storage. Designed to sit below the horizontal rail.



240907-(X)-BB

A - Amico Horizontal Rail
F - Fairfield Rail

E25 - DIN 25mm Rail

DW - Direct To Wall

V223 - 24" Surface Flatwall Spanning 2
Vertical Rails

Part#

* 240907-A-BB
* 240907-F-BB
* 240907-E25-BB

Part#

#* 240907-DW
^* 240907-V223-BB

^ 22.3" (56.6cm) between centers of two mounting points. Not compatible with Modular.

* Long lead time product.

Fasteners required to mount to wall are not provided.

Basket 28 x 5 x 3

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

28"^(71.1cm) x 5"^(12.7cm) x 3"^(7.6cm) (w x d x h)

Large basket designed to have room to accommodate items such as tubing, catheters, bandages etc. Designed to sit below the horizontal rail, or span across two vertical rails.



280503-(X)-BB

A - Amico Horizontal Rail
F - Fairfield Rail

E25 - DIN 25mm Rail

DW - Direct To Wall

V268 - 28" Amico Surface Flatwall
Spanning 2 Vertical Rails

Part#

* 280503-A-BB
* 280503-F-BB
* 280503-E25-BB

Part#

#* 280503-DW
^* 280503-V268-BB

* Long lead time product.

Fasteners required to mount to wall are not provided.

^ 26.81" (68.1cm) between centers of two mounting points. Not compatible with Modular.

Basket 28 x 5 x 6

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

28" ^(71.1cm) x 5" ^(12.7cm) x 6" ^(15.2cm) (w x d x h)

Storage basket designed to create maximum storage space for nurses and healthcare professionals.



280506-(X)-BB

A - Amico Horizontal Rail
F - Fairfield Rail

E25 - DIN 25mm Rail

DW - Direct To Wall

V268 - 28" Amico Surface Flatwall
Spanning 2 Vertical Rails

Part#

- * 280506-A-BB
- * 280506-F-BB
- * 280506-E25-BB

Part#

- #* 280506-DW
- ^* 280506-V268-BB

* Long lead time product.

Fasteners required to mount to wall are not provided.

^ 26.81" (68.1cm) between centers of two mounting points. Not compatible with Modular.

Basket 9 x 9 x 10 - Swivel Basket

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

9" ^(22.9cm) x 9" ^(22.9cm) x 10" ^(25.4cm) (w x d x h)

Medium size basket ideal to use as a patient or nurse trash bin when combined with a simple plastic disposable bag. Not compatible with Modular.



090910-(X)-BS

V - Amico Vertical Rail

Part#

090910-V-BS

Basket 10 x 6 x 9 - Swivel Basket

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

10" ^(25.4cm) x 6" ^(15.2cm) x 9" ^(22.9cm) (w x d x h)

Convenient storage basket with the addition of a swivel adapter that allows the basket to pivot side to side. Great for storage of gauze, oxygen masks, wipes, and other miscellaneous medical necessities. Not compatible with Modular.



100609-(X)-BS

V - Amico Vertical Rail

Part#

100609-V-BS

Basket 10 x 10 x 4 - Swivel Basket

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

10" (25.4cm) x 10" (25.4cm) x 4" (10.2cm) (w x d x h)

Versatile swivel basket which allows items to be within nurses reach while being kept organized. Not compatible with Modular.



101004-(X)-BS

V - Amico Vertical Rail

Part#

101004-V-BS

Basket 10 x 10 x 7 - Swivel Basket

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

10" (25.4cm) x 10" (25.4cm) x 7" (17.8cm) (w x d x h)

Medium size storage basket that helps organize and improve workflow. Not compatible with Modular.



101007-(X)-BS

V - Amico Vertical Rail

Part#

101007-V-BS

Basket 13 x 9 x 7 - Swivel Basket

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

13" (33.0cm) x 9" (22.9cm) x 7" (17.8cm) (w x d x h)

Large storage basket commonly used for patients belongings or other miscellaneous items, easily unhooks for transfer between rooms. Not compatible with Modular.



130907-(X)-BS

V - Amico Vertical Rail

Part#

130907-V-BS

Basket 18 x 9 x 7 - Raised Post Mount

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

18" (45.7cm) x 9" (22.9cm) x 7" (17.8cm) (w x d x h)

Large storage basket useful for those hard to store items that take up space on a shelf or in a cupboard. Designed to sit above the horizontal rail.



180907-(X)-BT

A - Amico Horizontal Rail

F - Fairfield Rail

E25 - DIN 25mm Rail

Part#

180907-A-BT

180907-F-BT

Part#

180907-E25-BT

Basket 19 x 5 x 3 - Raised Post Mount

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

19" (48.3cm) x 5" (12.7cm) x 3" (7.6cm) (w x d x h)

Medium storage basket accommodates larger items to free up space throughout the room. Designed to sit above the horizontal rail.



190503-(X)-BT

A - Amico Horizontal Rail

F - Fairfield Rail

E25 - DIN 25mm Rail

Part#

190503-A-BT

190503-F-BT

Part#

190503-E25-BT

Basket 24 x 5 x 3 - Raised Post Mount

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

24" (61.0cm) x 5" (12.7cm) x 3" (7.6cm) (w x d x h)

Large storage basket that can be moved and relocated easily on any part of the rail. Designed to sit above the horizontal rail.



240503-(X)-BT

A - Amico Horizontal Rail

F - Fairfield Rail

E25 - DIN 25mm Rail

Part#

240503-A-BT

240503-F-BT

Part#

240503-E25-BT

Basket 24 x 9 x 7 - Raised Post Mount

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

24" (61.0cm) x 9" (22.9cm) x 7" (17.8cm) (w x d x h)

Storage basket with ample room for larger items that require storage. Designed to sit above the horizontal rail.



240907-(X)-BT

A - Amico Horizontal Rail

F - Fairfield Rail

E25 - DIN 25mm Rail

Part#

* 240907-A-BT

* 240907-F-BT

Part#

* 240907-E25-BT

* Long lead time product.

Basket 28 x 5 x 3 - Raised Post Mount

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

28" (71.1cm) x 5" (12.7cm) x 3" (7.6cm) (w x d x h)

Large basket designed to accommodate items such as tubing, catheter, bandages etc. Designed to sit above the horizontal rail.



280503-(X)-BT

A - Amico Horizontal Rail

F - Fairfield Rail

E25 - DIN 25mm Rail

Part#

* 280503-A-BT

* 280503-F-BT

Part#

* 280503-E25-BT

* Long lead time product.

Basket 28 x 5 x 6 - Raised Post Mount

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

28" ^(71.1cm) x 5" ^(12.7cm) x 6" ^(15.2cm) (w x d x h)

Storage basket designed to sit above horizontal rails maximizing storage space for nurses and healthcare professionals.



280506-(X)-BT

A - Amico Horizontal Rail

F - Fairfield Rail

E25 - DIN 25mm Rail

Part#

- * 280506-A-BT
- * 280506-F-BT

Part#

- * 280506-E25-BT

* Long lead time product.

Basket Liner and Divider



DIV-101004-LNR



DIV-0503



DIV-097

Part#

Description

DIV-101004-LNR

Liner designed to fit a 10" x 10" x 4" (25.4cm x 25.4cm x 10.2cm) basket allowing separation of items. Ideal for smaller items.

DIV-0503

Basket item divider designed to fit a 5" (12.7cm) wide by 3" (7.6cm) deep basket for additional flexibility to fit your storage needs.

DIV-097

Basket item divider designed to fit a 9" (22.9 cm) wide by 7" (17.8 cm) deep basket for additional flexibility to fit your storage needs.

Folding Swivel Writing Shelf 6 x 6

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

6" ^(15.2cm) x 6" ^(15.2cm)

SHELFS-06FLD-(X)

Small work surface that allows pivoting and folding to fit in any space.



A - Amico Horizontal Rail

F - Fairfield Rail

V - Amico Vertical Rail

D - Ohmeda Dovetail Vertical Rail

G - Hill-Rom Vertical Rail

S - HSI Vertical Rail

E25 - DIN 25mm Rail

DW - Direct To Wall

M - Vertical Monitor Channel

Part#

- SHELFS-06FLD-A
- SHELFS-06FLD-F
- SHELFS-06FLD-V
- SHELFS-06FLD-D
- SHELFS-06FLD-G

Part#

- * SHELFS-06FLD-S
- SHELFS-06FLD-E25
- # SHELFS-06FLD-DW
- SHELFS-06FLD-M

* Long lead time product.

Fasteners required to mount to wall are not provided.

Folding Swivel Writing Shelf 12 x 12

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

12"^(30.5cm) X 12"^(30.5cm)

SHELFS-12FLD-(X)

Medium work surface that allows pivoting and folding to fit in any space.



A - Amico Horizontal Rail
F - Fairfield Rail
V - Amico Vertical Rail

D - Ohmeda Dovetail Vertical Rail
G - Hill-Rom Vertical Rail

S - HSI Vertical Rail
E25 - DIN 25mm Rail

DW - Direct To Wall
M - Vertical Monitor Channel

Part#

SHELFS-12FLD-A
SHELFS-12FLD-F
SHELFS-12FLD-V
SHELFS-12FLD-D
SHELFS-12FLD-G

Part#

* SHELFS-12FLD-S
SHELFS-12FLD-E25
#* SHELFS-12FLD-DW
* SHELFS-12FLD-M

Recommended maximum weight capacity: 5lbs (2.27kg)

* Long lead time product.

Fasteners required to mount to wall are not provided.

Folding Writing Shelf 20 x 12

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

19.5"^(49.5cm) X 12"^(30.5cm)

SHELFW-20FLD-(X)

Large work surface with smooth powder coated surfaces for easy cleaning.



A - Amico Horizontal Rail

F - Fairfield Rail

E25 - DIN 25mm Rail

V178 - 19" Amico Surface Flatwall
Spanning 2 Vertical Rails

Part#

SHELFW-20FLD-A
SHELFW-20FLD-F

Part#

SHELFW-20FLD-E25
^ SHELFW-20FLD-V178

Recommended maximum weight capacity: 20lbs (9.1kg)

^ 17.88"(45.4cm) between centers of two adapters.

Folding Writing Shelf 29 x 12

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

28.5"^(71.1cm) X 12"^(30.5cm)

SHELFW-29FLD-(X)

Extra large work surface with protective bumper for maximum durability.



A - Amico Horizontal Rail

F - Fairfield Rail

E25 - DIN 25mm Rail

V268 - 28" Amico Surface Flatwall
Spanning 2 Vertical Rails

Part#

SHELFW-29FLD-A
SHELFW-29FLD-F

Part#

SHELFW-29FLD-E25
^ SHELFW-29FLD-V268

^ 26.81"(68.1cm) between centers of two adapters.

Articulating Swivel Shelf

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

SHELFE-15-(X)

19.5" (49.5cm) X 12" (30.5cm)

Heavy duty shelf with 12"(30.5cm) extender arm for added flexibility (Weight rating: 15lbs).



V - Amico Vertical Rail

D - Ohmeda Dovetail Vertical Rail

MON - Vertical Monitor Channel

Part#

SHELFE-15-V

SHELFE-15-D

Part#

SHELFE-15-MON

Welch Allyn 767/777 Ophthalmoscope Mounting Plate

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

OPTHALMNT-(X)

Easy installation mounting plate for Welch Allyn 767 and 777 Oscope/Ophthalmoscope, available for virtually any rail system.



A - Amico Horizontal Rail

F - Fairfield Rail

V - Amico Vertical Rail

D - Ohmeda Dovetail Vertical Rail

G - Hill-Rom Vertical Rail

S - HSI Vertical Rail

E25 - DIN 25mm Rail

MON - Vertical Monitor Channel

Part#

OPTHALMNT-A

OPTHALMNT-F

OPTHALMNT-V

OPTHALMNT-D

OPTHALMNT-G

Part#

* OPTHALMNT-S

OPTHALMNT-E25

OPTHALMNT-MON

OPTHALMNT

* Long lead time product.

Welch Allyn Sure Temp/KleenSpec/Genius2 Tympanic Holder

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

SPECHLD-(X)

Designed to fit both Sure Temp (690, 692, 767, 777) thermometer, KleenSpec speculum dispenser, Genius 2 Tympanic or Filac 3000 thermometer.



A - Amico Horizontal Rail

F - Fairfield Rail

V - Amico Vertical Rail

D - Ohmeda Dovetail Vertical Rail

G - Hill-Rom Vertical Rail

S - HSI Vertical Rail

E25 - DIN 25mm Rail

R125 - 1.25" Diam. Pole Clamp

R200 - 2" Diam. Pole Clamp

MON - Vertical Monitor Channel

Part#

SPECHLD-A

SPECHLD-F

SPECHLD-V

SPECHLD-D

SPECHLD-G

Part#

* SPECHLD-S

SPECHLD-E25

SPECHLD-R125

SPECHLD-R200

SPECHLD-MON

* Long lead time product.

Exergen Instrument Holder Mount (134305 and 134306)

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

XG-INHLD-BKT-(X)

Designed to secure the Exergen 134305/134306 keyless locking wall mount to various Amico wall rails. Customer to supply the Exergen 134305/134306 keyless locking wall mount.



A - Amico Horizontal Rail F - Fairfield Rail V - Amico Vertical Rail

Part#

- XG-INHLD-BKT-A
- XG-INHLD-BKT-F
- XG-INHLD-BKT-V

BP Cuff Basket

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

AM-CUFBSK-(X)

Basket designed to accommodate all cuff and inflation systems.



A - Amico Horizontal Rail D - Ohmeda Dovetail Vertical Rail E25 - DIN 25mm Rail R200 - 2" Diam. Pole Clamp
F - Fairfield Rail G - Hill-Rom Vertical Rail R125 - 1,25" Diam. Pole Clamp M - Vertical Monitor Channel
V - Amico Vertical Rail S - HSI Vertical Rail

Part#

- AM-CUFBSK-A
- AM-CUFBSK-F
- AM-CUFBSK-V
- AM-CUFBSK-D
- AM-CUFBSK-G
- * AM-CUFBSK-S

Part#

- AM-CUFBSK-E25
- AM-CUFBSK-R125
- AM-CUFBSK-R200
- AM-CUFBSK-M
- # AM-CUFBSK

Fasteners required to mount to wall are not provided.

* Long lead time product.

Amico Aneroid Manometer With Cuff Basket

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

AM-(X)L-LF2118

Amico Aneroid Manometer including cuff basket and the ability to swivel for optimal viewing. System includes gauge, inflation system, adult cuff, and cuff basket.



A - Amico Horizontal Rail V - Amico Vertical Rail S - HSI Vertical Rail M - Vertical Monitor Channel
F - Fairfield Rail G - Hill-Rom Vertical Rail E25 - DIN 25mm Rail

Part#

- AM-AL-LF2118
- AM-FL-LF2118
- AM-VL-LF2118
- AM-GL-LF2118

Part#

- * AM-SL-LF2118
- AM-EL-LF2118
- AM-ML-LF2118
- AM-XL-LF2118

* Long lead time product.

Welch Allyn 767 Aneroid Mounting Plate

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

AM-767MNT-(X)

Compatible with all Welch Allyn 767 Aneroid Manometers to mount to various styles of accessory rails.



A - Amico Horizontal Rail	D - Ohmeda Dovetail Vertical Rail	E25 - DIN 25mm Rail	R200 - 2" Diam. Pole Clamp
F - Fairfield Rail	G - Hill-Rom Vertical Rail	R125 - 1.25" Diam. Pole Clamp	MON - Vertical Monitor Channel
V - Amico Vertical Rail	S - HSI Vertical Rail		

Part#

AM-767MNT-A
AM-767MNT-F
AM-767MNT-V
AM-767MNT-D
AM-767MNT-G

Part#

* AM-767MNT-S
AM-767MNT-E25
AM-767MNT-R125
AM-767MNT-R200
AM-767MNT-MON

* Long lead time product.

Welch Allyn 5091 Aneroid Mounting Plate

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

AM-5091MNT-(X)

Mounting plate that works with the Welch Allyn 5091 Aneroid Manometer to mount to various styles of accessory rails.



A - Amico Horizontal Rail	D - Ohmeda Dovetail Vertical Rail	E25 - DIN 25mm Rail	R200 - 2" Diam. Pole Clamp
F - Fairfield Rail	G - Hill-Rom Vertical Rail	R125 - 1.25" Diam. Pole Clamp	MON - Vertical Monitor Channel
V - Amico Vertical Rail	S - HSI Vertical Rail		

Part#

AM-5091MNT-A
AM-5091MNT-F
AM-5091MNT-V
AM-5091MNT-D
AM-5091MNT-G

Part#

* AM-5091MNT-S
AM-5091MNT-E25
AM-5091MNT-R125
AM-5091MNT-R200
AM-5091MNT-MON

* Long lead time product.

TrimLine AM4300 Aneroid Mounting Plate

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

AM-4300MNT-(X)

Mounting plate compatible with the TrimLine AM4300 Aneroid Manometer to mount device on various styles of accessory rails.



A - Amico Horizontal Rail	D - Ohmeda Dovetail Vertical Rail	E25 - DIN 25mm Rail	R200 - 2" Diam. Pole Clamp
F - Fairfield Rail	G - Hill-Rom Vertical Rail	R125 - 1.25" Diam. Pole Clamp	MON - Vertical Monitor Channel
V - Amico Vertical Rail	S - HSI Vertical Rail		

Part#

AM-4300MNT-A
AM-4300MNT-F
AM-4300MNT-V
AM-4300MNT-D
AM-4300MNT-G

Part#

* AM-4300MNT-S
AM-4300MNT-E25
AM-4300MNT-R125
AM-4300MNT-R200
AM-4300MNT-MON

* Long lead time product.

WA Baum or ADC 750W Aneroid Mounting Plate

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

AM-750WMNT-(X)

WA Baum or ADC 750W Aneroid mounting plate with adapter. Allows equipment to be easily moved to maximize workflow.



A - Amico Horizontal Rail	D - Ohmeda Dovetail Vertical Rail	E25 - DIN 25mm Rail	R200 - 2" Diam. Pole Clamp
F - Fairfield Rail	G - Hill-Rom Vertical Rail	R125 - 1.25" Diam. Pole Clamp	MON - Vertical Monitor Channel
V - Amico Vertical Rail	S - HSI Vertical Rail		

Part#

AM-750WMNT-A
AM-750WMNT-F
AM-750WMNT-V
AM-750WMNT-D
AM-750WMNT-G

Part#

* AM-750WMNT-S
AM-750WMNT-E25
AM-750WMNT-R125
AM-750WMNT-R200
AM-750WMNT-MON

* Long lead time product.

Vacuum Bottle Slide

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

SLIDE-(X)

Universal suction canister ring holder. Designed to be easily moved and easily installed.



A - Amico Horizontal Rail	D - Ohmeda Dovetail Vertical Rail	E25 - DIN 25mm Rail	MON - Vertical Monitor Channel
F - Fairfield Rail	G - Hill-Rom Vertical Rail	R125 - 1.25" Diam. Pole Clamp	
V - Amico Vertical Rail	S - HSI Vertical Rail	R200 - 2" Diam. Pole Clamp	

Part#

SLIDE-A
SLIDE-F
SLIDE-V
SLIDE-D
SLIDE-G

* SLIDE-S

Part#

SLIDE-E25
SLIDE-R125
SLIDE-R200
SLIDE-MON
SLIDE

* Long lead time product.

Plastic Suction Canister Holder 1200cc/mL

VR-CNHL-1200

Holder for the 1200cc disposable suction canister
(Slide not included).



Part#

VR-CNHL-1200

Plastic Suction Canister Holder 2000cc/mL

VR-CNHL-2000

Holder for the 2000cc & 3000cc disposable suction
canister (Slide not included).



Part#

VR-CNHL-2000

Plastic Suction Canister 1200cc/mL

VR-CAN-1200

1200cc disposable suction canister, ready to use out
of the box.



Part#

VR-CAN-1200

Plastic Suction Canister 2000cc/mL

VR-CAN-2000

2000cc disposable suction canister, ready to use out
of the box.



Part#

VR-CAN-2000

Plastic Suction Canister 3000cc/mL

VR-CAN-3000

3000cc disposable suction canister, ready to use out
of the box.



Part#

* VR-CAN-3000

* Long lead time product.

Blender Bracket

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

BLBRKT-(X)

Anodized aluminum prevents scratches and slows down wear, compatible with all makes of O₂ Blenders. Mounting options available with any rail or power column.



A - Amico Horizontal Rail	D - Ohmeda Dovetail Vertical Rail	E25 - DIN 25mm Rail	MON - Vertical Monitor Channel
F - Fairfield Rail	G - Hill-Rom Vertical Rail	R125 - 1.25" Diam. Pole Clamp	
V - Amico Vertical Rail	S - HSI Vertical Rail	R200 - 2" Diam. Pole Clamp	

Part#

BLBRKT-A
BLBRKT-F
BLBRKT-V
BLBRKT-D
BLBRKT-G
* BLBRKT-S

Part#

BLBRKT-E25
BLBRKT-R125
BLBRKT-R200
BLBRKT-MON
BLBRKT
BLBRKT-PC1

* Long lead time product.

Humidifier Bracket

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

HUMBRKT-(X)

Fisher & Paykel heated humidifier mount with bottom support for secure mounting.



A - Amico Horizontal Rail	D - Ohmeda Dovetail Vertical Rail	E25 - DIN 25mm Rail	MON - Vertical Monitor Channel
F - Fairfield Rail	G - Hill-Rom Vertical Rail	R125 - 1.25" Diam. Pole Clamp	
V - Amico Vertical Rail	S - HSI Vertical Rail	R200 - 2" Diam. Pole Clamp	

Part#

HUMBRKT-A
HUMBRKT-F
HUMBRKT-V
HUMBRKT-D
HUMBRKT-G
* HUMBRKT-S

Part#

HUMBRKT-E25
HUMBRKT-R125
HUMBRKT-R200
HUMBRKT-MON
HUMBRKT

* Long lead time product.

Supply Clip - Single

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

SUPCLIP-1-(X)

Ideal for holding single use packaged items with quick pull release. Not compatible with Modular.



F - Fairfield Rail	G - Hill-Rom Vertical Rail	E25 - DIN 25mm Rail	R200 - 2" Diam. Pole Clamp
V - Amico Vertical Rail	S - HSI Vertical Rail	R125 - 1.25" Diam. Pole Clamp	MON - Vertical Monitor Channel
D - Ohmeda Dovetail Vertical Rail			

Part#

SUPCLIP-1-F
SUPCLIP-1-V
SUPCLIP-1-D
SUPCLIP-1-G
* SUPCLIP-1-S

Part#

SUPCLIP-1-E25
SUPCLIP-1-R125
SUPCLIP-1-R200
SUPCLIP-1-MON

* Long lead time product.

Supply Clip - Dual

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

SUPCLIP-2-(X)

Ideal for holding single use packaged items with quick pull release - Dual clips included for two in one organization.



A - Amico Horizontal Rail
F - Fairfield Rail
D - Ohmeda Dovetail Vertical Rail

G - Hill-Rom Vertical Rail
S - HSI Vertical Rail

E25 - DIN 25mm Rail
R125 - 1.25" Diam. Pole Clamp

R200 - 2" Diam. Pole Clamp
MON - Vertical Monitor Channel

Part#

SUPCLIP-2-A

- * SUPCLIP-2-F
- * SUPCLIP-2-D
- * SUPCLIP-2-G
- * SUPCLIP-2-S

Part#

- * SUPCLIP-2-E25
- * SUPCLIP-2-R125
- * SUPCLIP-2-R200
- * SUPCLIP-2-MON

* Long lead time product.

Tubing Clip - Single

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

TUBECLIP-1-(X)

Useful for hoses, electrical , and nurse call cords, also neatly organizes O₂ Cannulas and tubes. Not compatible with Modular.



F - Fairfield Rail
V - Amico Vertical Rail
D - Ohmeda Dovetail Vertical Rail

G - Hill-Rom Vertical Rail
S - HSI Vertical Rail

E25 - DIN 25mm Rail
R125 - 1.25" Diam. Pole Clamp

R200 - 2" Diam. Pole Clamp
MON - Vertical Monitor Channel

Part#

TUBECLIP-1-F
TUBECLIP-1-V
TUBECLIP-1-D
TUBECLIP-1-G

- * TUBECLIP-1-S

Part#

TUBECLIP-1-E25
TUBECLIP-1-R125
TUBECLIP-1-R200
TUBECLIP-1-MON

* Long lead time product.

Tubing Clip - Double

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

TUBECLIP-2-(X)

Useful for hoses, electrical , and nurse call cords, also neatly organizes O₂ Cannulas and tubes - includes dual clips.



A - Amico Horizontal Rail
F - Fairfield Rail
D - Ohmeda Dovetail Vertical Rail

G - Hill-Rom Vertical Rail
S - HSI Vertical Rail

E25 - DIN 25mm Rail
R125 - 1.25" Diam. Pole Clamp

R200 - 2" Diam. Pole Clamp
MON - Vertical Monitor Channel

Part#

TUBECLIP-2-A
* TUBECLIP-2-F
* TUBECLIP-2-D
* TUBECLIP-2-G
* TUBECLIP-2-S

Part#

- * TUBECLIP-2-E25
- * TUBECLIP-2-R125
- * TUBECLIP-2-R200
- * TUBECLIP-2-MON

* Long lead time product.

Utility Clip

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

UTILCLIP-(X)

Versatile solution in rooms for storing tubing, BP cuff, cables, etc.



A - Amico Horizontal Rail	D - Ohmeda Dovetail Vertical Rail	E25 - DIN 25mm Rail	MON - Vertical Monitor Channel
F - Fairfield Rail	G - Hill-Rom Vertical Rail	R125 - 1.25" Diam. Pole Clamp	
V - Amico Vertical Rail	S - HSI Vertical Rail	R200 - 2" Diam. Pole Clamp	

Part#

UTILCLIP-A
UTILCLIP-F
UTILCLIP-V
UTILCLIP-D
UTILCLIP-G
* UTILCLIP-S

Part#

UTILCLIP-E25
UTILCLIP-R125
UTILCLIP-R200
UTILCLIP-MON
UTILCLIP

* Long lead time product.

Cord Wrap

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

CORDWRAP-(X)

Ideal for managing wires and cords to create a cleaner working space.



A - Amico Horizontal Rail	D - Ohmeda Dovetail Vertical Rail	E25 - DIN 25mm Rail	MON - Vertical Monitor Channel
F - Fairfield Rail	G - Hill-Rom Vertical Rail	R125 - 1.25" Diam. Pole Clamp	
V - Amico Vertical Rail	S - HSI Vertical Rail	R200 - 2" Diam. Pole Clamp	

Part#

CORDWRAP-A
CORDWRAP-F
CORDWRAP-V
CORDWRAP-D
CORDWRAP-G
* CORDWRAP-S

Part#

CORDWRAP-E25
CORDWRAP-R125
CORDWRAP-R200
CORDWRAP-MON
CORDWRAP

* Long lead time product.

Chart Holder

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

CHARTHLD-(X)

Efficient patient chart organizer. Helps keep room looking neat and well organized.



A - Amico Horizontal Rail	D - Ohmeda Dovetail Vertical Rail	E25 - DIN 25mm Rail	MON - Vertical Monitor Channel
F - Fairfield Rail	G - Hill-Rom Vertical Rail	R125 - 1.25" Diam. Pole Clamp	
V - Amico Vertical Rail	S - HSI Vertical Rail	R200 - 2" Diam. Pole Clamp	

Part#

CHARTHLD-A
CHARTHLD-F
CHARTHLD-V
CHARTHLD-D
CHARTHLD-G
* CHARTHLD-S

Part#

CHARTHLD-E25
CHARTHLD-R125
CHARTHLD-R200
CHARTHLD-MON
CHARTHLD

* Long lead time product.

Resuscitation Bag Holder

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

RESUSHLD-(X)

Dual hook design provides an ideal space to hold
resuscitation bag.



A - Amico Horizontal Rail
F - Fairfield Rail

V - Amico Vertical Rail
D - Ohmeda Dovetail Vertical Rail

G - Hill-Rom Vertical Rail
S - HSI Vertical Rail

E25 - DIN 25mm Rail
MON - Vertical Monitor Channel

Part#

RESUSHLD-A
RESUSHLD-F
RESUSHLD-V
RESUSHLD-D
RESUSHLD-G

Part#

* RESUSHLD-S
RESUSHLD-E25
RESUSHLD-MON
RESUSHLD

* Long lead time product.

E-Cylinder Tank Holder

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

ECYLHLDR-(X)

E-cylinder holder with secure locking mechanism,
available in a white powder coat finish.



F - Fairfield Rail

E25 - DIN 25mm Rail

Part#

ECYLHLDR-F

Part#

ECYLHLDR-E25

Fit type E cylinder only.

Daniels Sharpsmart Holder

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

DA-SHARPSHLD-(X)

Designed to fit Daniels Sharpsmart Sharps Holder
and can be mounted to any rail system. Designed to
work with Daniels product code 4005007000.



A - Amico Horizontal Rail
F - Fairfield Rail
V - Amico Vertical Rail

D - Ohmeda Dovetail Vertical Rail
G - Hill-Rom Vertical Rail
S - HSI Vertical Rail

E25 - DIN 25mm Rail
R125 - 1.25" Diam. Pole Clamp

R200 - 2" Diam. Pole Clamp
MON - Vertical Monitor Channel

Part#

* DA-SHARPSHLD-A
* DA-SHARPSHLD-F
* DA-SHARPSHLD-V
* DA-SHARPSHLD-D
* DA-SHARPSHLD-G

Part#

* DA-SHARPSHLD-S
* DA-SHARPSHLD-E25
* DA-SHARPSHLD-R125
* DA-SHARPSHLD-R200
* DA-SHARPSHLD-MON

* Long lead time product.

Universal Sharps Holder

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

SHARPSHLD-(X)

Designed to be easily mounted and ready to use out of the box. Created to work with any brand of sharps container with footprint max. at 10"x10".



A - Amico Horizontal Rail	D - Ohmeda Dovetail Vertical Rail	E25 - DIN 25mm Rail	MON - Vertical Monitor Channel
F - Fairfield Rail	G - Hill-Rom Vertical Rail	R125 - 1.25" Diam. Pole Clamp	
V - Amico Vertical Rail	S - HSI Vertical Rail	R200 - 2" Diam. Pole Clamp	

Part#

SHARPSHLD-A
SHARPSHLD-F
SHARPSHLD-V
SHARPSHLD-D
SHARPSHLD-G
* SHARPSHLD-S

Part#

SHARPSHLD-E25
SHARPSHLD-R125
SHARPSHLD-R200
SHARPSHLD-MON
SHARPSHLD

* Long lead time product.

Glove Box Holder - Single

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

GLBOXHLD-1-(X)

Versatile storage for glove boxes or tissue boxes, can be mounted onto various style rails or direct to wall.



A - Amico Horizontal Rail	D - Ohmeda Dovetail Vertical Rail	E25 - DIN 25mm Rail	MON - Vertical Monitor Channel
F - Fairfield Rail	G - Hill-Rom Vertical Rail	R125 - 1.25" Diam. Pole Clamp	
V - Amico Vertical Rail	S - HSI Vertical Rail	R200 - 2" Diam. Pole Clamp	

Part#

GLBOXHLD-1-A
GLBOXHLD-1-F
GLBOXHLD-1-V
GLBOXHLD-1-D
GLBOXHLD-1-G
* GLBOXHLD-1-S

Part#

GLBOXHLD-1-E25
GLBOXHLD-1-R125
GLBOXHLD-1-R200
GLBOXHLD-1-MON
GLBOXHLD-1

* Long lead time product.

Quad (4 Hook) IV Hook

IVHOOK4

Quad IV Hook with 9"(22.9cm) of height adjustment allowing users to maximize small spaces.



Part#

IVHOOK4

Required a clearance of 5"(12.7cm) from wall. Does not include mounting adapter, see page 34.

Duplex (2 Hook) IV Hook

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

IVHOOK2-(X)

Dual IV Hook with 9"(22.9cm) of height adjustment
designed to keep IV bags accessible.



A - Amico Horizontal Rail	V - Amico Vertical Rail	G - Hill-Rom Vertical Rail	E25 - DIN 25mm Rail
F - Fairfield Rail	D - Ohmeda Dovetail Vertical Rail	S - HSI Vertical Rail	MON - Vertical Monitor Channel

Part#

IVHOOK2-A
IVHOOK2-F
IVHOOK2-V
IVHOOK2-D
IVHOOK2-G

Part#

* IVHOOK2-S
IVHOOK2-E25
IVHOOK2-MON
IVHOOK2

* Long lead time product.

36" Duplex (2 Hook) IV Hook

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

IVHOOK2-36-(X)

36"(91.4cm) dual IV Hook with dual mounting points
for added stability provides a maximum flexibility on
placement for IV bags.



A75 - Amico (7.5" c-c) Spans 2 Horizontal Rails	A18 - (18" c-c) Spans 2 Amico Horizontal Rails 18" Apart	V - Amico Vertical Rail
---	--	-------------------------

Part#

* IVHOOK2-36-A75
* IVHOOK2-36-A18

Part#

* IVHOOK2-36
* IVHOOK2-36-V

* Long lead time product.

Infusion Pump Mount

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

INFPMP-(X)

Durable infusion pump mount with 1" diameter
swivel aluminum pole for pumps and a quad IV bag
holder.



V - Amico Vertical Rail	S - HSI Vertical Rail	A18 - (18" c-c) Spans 2 Amico Horizontal Rails 18" Apart	A75 - Amico (7.5" c-c) Spans 2 Horizontal Rails
F - Fairfield Rail	E25 - DIN 25mm Rail		
G - Hill-Rom Vertical Rail			

Part#

INFPMP-V
INFPMP-F
INFPMP-G
* INFPMP-S

Part#

INFPMP-E25
† INFPMP-A18
INFPMP-A75
INFPMP

23" mounting surface (23.75" between adapters c-c) on all models except INFPMP-A18 with surface at 17" (17.75" between adapters c-c) and INFPMP-A75 with surface at 6.5" (7.25" between adapters c-c).

* Long lead time product.






† Equivalent to previous part# INFPMP-A.


Infusion Extension

Replace (X) with selected extension type, (Y) with hook style and (Z) with selected adapter style. Shipped fully assembled with selected extensions, hooks and adapter style.

Infusion Extension with 12" up post and 12" down post for secure mounting of multiple IV bags and infusion pumps.

INE-09(XX)-12-12-(YY)-(Z)

Extension (XX)	IV Hooks (YY)	Adapter Style (Z)
 <p>00 9" Extension</p>  <p>09 9"x9" Folding Extension</p>	<p>00 No Hooks</p> <p>F4 4 Hooks</p>	 <p>AT Amico Horizontal Mount-Top</p>  <p>V Amico Vertical Mount</p>  <p>M Monitor Channel Mount</p>



9" Extension by Adapter Style

Across 2 Amico Horizontal Rails (Top Mount) - AT

Part#

INE-0900-12-12-00-AT
INE-0900-12-12-F4-AT

Amico Vertical Rail - V

Part#

INE-0900-12-12-00-V
INE-0900-12-12-F4-V

Vertical Monitor Channel - M

Part#

INE-0900-12-12-00-M
INE-0900-12-12-F4-M

9x9" Extension by Adapter Style

Across 2 Amico Horizontal Rails (Top Mount) - AT

Part#

INE-0909-12-12-00-AT
INE-0909-12-12-F4-AT

Amico Vertical Rail - V

Part#

INE-0909-12-12-00-V
INE-0909-12-12-F4-V

Vertical Monitor Channel - M

Part#

INE-0909-12-12-00-M
INE-0909-12-12-F4-M

Standard in white for the arm only. Other color options available with additional cost. Please call for details.

Equipment Post Straight 1" x 10"

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

POST-S-1010-(X)

1" (2.54cm) x 10" (25.4cm)

Designed for lightweight pole-mounted device, like syringe pump and transducers, to stay accessible and tucked away.



A - Amico Horizontal Rail	D - Ohmeda Dovetail Vertical Rail	E25 - DIN 25mm Rail	MON - Vertical Monitor Channel
F - Fairfield Rail	G - Hill-Rom Vertical Rail	DW - Direct To Wall	
V - Amico Vertical Rail	S - HSI Vertical Rail		

Part#

POST-S-1010-A
POST-S-1010-F
POST-S-1010-V
POST-S-1010-D
POST-S-1010-G

Part#

* POST-S-1010-S
POST-S-1010-E25
POST-S-1010-DW
POST-S-1010-MON
POST-S-1010

* Long lead time product.

Fasteners required to mount to wall are not provided.

Equipment Post Offset 1" x 10"

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

POST-O-1010-(X)

1" (2.54cm) x 10" (25.4cm)

Designed for lightweight pole-mounted device, like syringe pump and transducers, offset design allows for easy access to plugs and clamps.



A - Amico Horizontal Rail	D - Ohmeda Dovetail Vertical Rail	E25 - DIN 25mm Rail	MON - Vertical Monitor Channel
F - Fairfield Rail	G - Hill-Rom Vertical Rail	DW - Direct To Wall	
V - Amico Vertical Rail	S - HSI Vertical Rail		

Part#

POST-O-1010-A
POST-O-1010-F
POST-O-1010-V
POST-O-1010-D
POST-O-1010-G

Part#

* POST-O-1010-S
POST-O-1010-E25
POST-O-1010-DW
POST-O-1010-MON
POST-O-1010

* Long lead time product.

Fasteners required to mount to wall are not provided.

Equipment Post Offset 1/2" x 3"

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

POST-O-0503-(X)

1/2" (1.27cm) x 3" (7.62cm)

Perfect for any lightweight 1/2" pole-mounting device or as a hook for patient particulars or supplies.



A - Amico Horizontal Rail	D - Ohmeda Dovetail Vertical Rail	E25 - DIN 25mm Rail	MON - Vertical Monitor Channel
F - Fairfield Rail	G - Hill-Rom Vertical Rail	DW - Direct To Wall	
V - Amico Vertical Rail	S - HSI Vertical Rail		

Part#

POST-O-0503-A
POST-O-0503-F
POST-O-0503-V
POST-O-0503-D
POST-O-0503-G

Part#

* POST-O-0503-S
POST-O-0503-E25
POST-O-0503-DW
POST-O-0503-MON
POST-O-0503

* Long lead time product.

Fasteners required to mount to wall are not provided.

Universal Barcode Scanner Holder

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

BARCD-(X)

Universal bar code scanner holder designed to accommodate most tethered scanners.



A - Amico Horizontal Rail	D - Ohmeda Dovetail Vertical Rail	S - HSI Vertical Rail	V - Amico Vertical Rail
F - Fairfield Rail	G - Hill-Rom Vertical Rail	MON - Vertical Monitor Channel	

Part#

BARCD-A
BARCD-F
BARCD-V
BARCD-D

Part#

BARCD-G
* BARCD-S
BARCD-MON

* Long lead time product.

Power Supply Mount

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

PWRHLD-(X)

Provides secure mounting for most power supplies
for laptops, patient monitors and LCD monitors.



A - Amico Horizontal Rail
F - Fairfield Rail
V - Amico Vertical Rail
D - Ohmeda Dovetail Vertical Rail

G - Hill-Rom Vertical Rail
S - HSI Vertical Rail
E25 - DIN 25mm Rail
R125 - 1.25" Diam. Pole Clamp

R200 - 2" Diam. Pole Clamp
MON - Vertical Monitor Channel
PLV - Adj. Pole Clamp for vert. pole w/
0.5-1.5" diam.

PLH - Adj. Pole Clamp for horiz. pole at
0.5-1.5" diam./DIN rail

Part#

PWRHLD-A
PWRHLD-F
PWRHLD-V
PWRHLD-D
PWRHLD-G
* PWRHLD-S
PWRHLD-E25

Part#

PWRHLD-R125
PWRHLD-R200
PWRHLD-MON
PWRHLD-PLV
PWRHLD-PLH
PWRHLD

* Long lead time product.

Power Supply Mount - Dual

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

PWRHLD2-(X)

Holds two power supplies at a time for most laptops,
patient monitors and LCD monitors.



A - Amico Horizontal Rail
F - Fairfield Rail
V - Amico Vertical Rail

D - Ohmeda Dovetail Vertical Rail
G - Hill-Rom Vertical Rail
S - HSI Vertical Rail

E25 - DIN 25mm Rail
R125 - 1.25" Diam. Pole Clamp
R200 - 2" Diam. Pole Clamp

MON - Vertical Monitor Channel

Part#

PWRHLD2-A
PWRHLD2-F
PWRHLD2-V
PWRHLD2-D
PWRHLD2-G
* PWRHLD2-S

Part#

PWRHLD2-E25
PWRHLD2-R125
PWRHLD2-R200
PWRHLD2-MON
PWRHLD2

* Long lead time product.

Single Gas Blocks - Oxygen

Replace (X) with selected adapter and (Y) with selected outlet. Shipped fully assembled with selected adapter style

Single Gas Blocks for Oxygen to add flexibility in moving your gas distribution closer to the patient, compatible with all gas outlets.



G(X)S-PL(Y)PL-DMUO

Adapter Style (X)			Outlet Type (Y)		
A Amico Horizontal	F Fairfield	V Vertical	DC DISS	OF Ohmeda Compatible	CF Chemetron Compatible
E1 DIN 25mm	W Wall / Slide		PF Puritan Bennett Compatible	MF Medstar Compatible	F4 1/4" FNPT Compatible

Rendering Photo Not Available

By Adapter Style

Amico Horizontal Rail - A

Part#

GAS-PLDCPL-DMUO
GAS-PLOFPL-DMUO
GAS-PLCFPL-DMUO
GAS-PLPFPL-DMUO
GAS-PLMFPL-DMUO
GAS-PLF4PL-DMUO

Amico Vertical Rail - V

Part#

GVS-PLDCPL-DMUO
GVS-PLOFPL-DMUO
GVS-PLCFPL-DMUO
GVS-PLPFPL-DMUO
GVS-PLMFPL-DMUO
GVS-PLF4PL-DMUO

Direct To Wall - W

Part#

GWS-PLDCPL-DMUO
GWS-PLOFPL-DMUO
GWS-PLCFPL-DMUO
GWS-PLPFPL-DMUO
GWS-PLMFPL-DMUO
GWS-PLF4PL-DMUO

Fairfield Rail - F

Part#

‡ GFS-PLDCPL-DMUO
‡ GFS-PLOFPL-DMUO
‡ GFS-PLCFPL-DMUO
‡ GFS-PLPFPL-DMUO
‡ GFS-PLMFPL-DMUO
‡ GFS-PLF4PL-DMUO

DIN 25mm Rail - E1

Part#

GE1S-PLDCPL-DMUO
GE1S-PLOFPL-DMUO
GE1S-PLCFPL-DMUO
GE1S-PLPFPL-DMUO
GE1S-PLMFPL-DMUO
GE1S-PLF4PL-DMUO

Built to order. Call for lead time.

Add \$7 to substitute a PL (plug) for PT (90° Elbow Down DISS with check).

Substitute U with I for ISO color codes at no extra charge.

‡ Add \$30 and change 'GFS' to 'GFPS' to substitute Fairfield CAM adapter with push button style.

Single Gas Blocks - Medical Air

Replace (X) with selected adapter and (Y) with selected outlet. Shipped fully assembled with selected adapter style

Single Gas Blocks for Medical Air to add flexibility in moving your gas distribution closer to the patient, compatible with all gas outlets.



G(X)S-PL(Y)PL-DMUA

Adapter Style (X)			Outlet Type (Y)		
A Amico Horizontal	F Fairfield	V Vertical	DC DISS	OF Ohmeda Compatible	CF Chemetron Compatible
E1 DIN 25mm	W Wall / Slide		PF Puritan Bennett Compatible	MF Medstar Compatible	F4 1/4" FNPT Compatible

By Adapter Style

Amico Horizontal Rail - A

Part#

GAS-PLDCPL-DMUA
GAS-PLOFPL-DMUA
GAS-PLCFPL-DMUA
GAS-PLPFPL-DMUA
GAS-PLMFPL-DMUA
GAS-PLF4PL-DMUA

Amico Vertical Rail - V

Part#

GVS-PLDCPL-DMUA
GVS-PLOFPL-DMUA
GVS-PLCFPL-DMUA
GVS-PLPFPL-DMUA
GVS-PLMFPL-DMUA
GVS-PLF4PL-DMUA

Direct To Wall - W

Part#

GWS-PLDCPL-DMUA
GWS-PLOFPL-DMUA
GWS-PLCFPL-DMUA
GWS-PLPFPL-DMUA
GWS-PLMFPL-DMUA
GWS-PLF4PL-DMUA

Fairfield Rail - F

Part#

‡ GFS-PLDCPL-DMUA
‡ GFS-PLOFPL-DMUA
‡ GFS-PLCFPL-DMUA
‡ GFS-PLPFPL-DMUA
‡ GFS-PLMFPL-DMUA
‡ GFS-PLF4PL-DMUA

DIN 25mm Rail - E1

Part#

GE1S-PLDCPL-DMUA
GE1S-PLOFPL-DMUA
GE1S-PLCFPL-DMUA
GE1S-PLPFPL-DMUA
GE1S-PLMFPL-DMUA
GE1S-PLF4PL-DMUA

Built to order. Call for lead time.

Add \$7 to substitute a PL (plug) for PT (90° Elbow Down DISS with check).


Substitute U with I for ISO color codes at no extra charge.

‡ Add \$30 and change 'GFS' to 'GFPS' to substitute Fairfield CAM adapter with push button style.

Single Gas Blocks - Vacuum

Replace (X) with selected adapter and (Y) with selected outlet. Shipped fully assembled with selected adapter style

Single Gas Blocks for vacuum to add flexibility in moving your gas distribution closer to the patient, compatible with all gas outlets.



By Adapter Style

Amico Horizontal Rail - A

Part#

GAS-PLDCPL-DMUV
GAS-PLOFPL-DMUV
GAS-PLCFPL-DMUV
GAS-PLPFPL-DMUV
GAS-PLMFPL-DMUV
GAS-PLF4PL-DMUV

Amico Vertical Rail - V

Part#

GVS-PLDCPL-DMUV
GVS-PLOFPL-DMUV
GVS-PLCFPL-DMUV
GVS-PLPFPL-DMUV
GVS-PLMFPL-DMUV
GVS-PLF4PL-DMUV

Direct To Wall - W

Part#

GWS-PLDCPL-DMUV
GWS-PLOFPL-DMUV
GWS-PLCFPL-DMUV
GWS-PLPFPL-DMUV
GWS-PLMFPL-DMUV
GWS-PLF4PL-DMUV

Fairfield Rail - F

Part#

‡ GFS-PLDCPL-DMUV
‡ GFS-PLOFPL-DMUV
‡ GFS-PLCFPL-DMUV
‡ GFS-PLPFPL-DMUV
‡ GFS-PLMFPL-DMUV
‡ GFS-PLF4PL-DMUV

DIN 25mm Rail - E1

Part#

GE1S-PLDCPL-DMUV
GE1S-PLOFPL-DMUV
GE1S-PLCFPL-DMUV
GE1S-PLPFPL-DMUV
GE1S-PLMFPL-DMUV
GE1S-PLF4PL-DMUV

Built to order. Call for lead time.

Add \$7 to substitute a PL (plug) for PT (90° Elbow Down DISS with check).

Substitute U with I for ISO color codes at no extra charge.

‡ Add \$30 and change 'GFS' to 'GFPS' to substitute Fairfield CAM adapter with push button style.

Double Gas Blocks - Oxygen

Replace (X) with selected adapter and (Y) with selected outlet. Shipped fully assembled with selected adapter style

Cost effective alternative to adding new oxygen outlets or moving outlets, suited for all gases.



G(X)D-PL(Y)PL-DMUO

Adapter Style (X)			Outlet Type (Y)		
A Amico Horizontal	F Fairfield	V Vertical	DC DISS	OF Ohmeda Compatible	CF Chemetron Compatible
E1 DIN 25mm	W Wall / Slide		PF Puritan Bennett Compatible	MF Medstar Compatible	F4 1/4" FNPT Compatible

By Adapter Style

Amico Horizontal Rail - A

Part#

GAD-PLDCPL-DMUO
GAD-PLOFPL-DMUO
GAD-PLCFPL-DMUO
GAD-PLPFPL-DMUO
GAD-PLMFPL-DMUO
GAD-PLF4PL-DMUO

Amico Vertical Rail - V

Part#

GVD-PLDCPL-DMUO
GVD-PLOFPL-DMUO
GVD-PLCFPL-DMUO
GVD-PLPFPL-DMUO
GVD-PLMFPL-DMUO
GVD-PLF4PL-DMUO

Direct To Wall - W

Part#

GWD-PLDCPL-DMUO
GWD-PLOFPL-DMUO
GWD-PLCFPL-DMUO
GWD-PLPFPL-DMUO
GWD-PLMFPL-DMUO
GWD-PLF4PL-DMUO

Fairfield Rail - F

Part#

‡ GFD-PLDCPL-DMUO
‡ GFD-PLOFPL-DMUO
‡ GFD-PLCFPL-DMUO
‡ GFD-PLPFPL-DMUO
‡ GFD-PLMFPL-DMUO
‡ GFD-PLF4PL-DMUO

DIN 25mm Rail - E1

Part#

GE1D-PLDCPL-DMUO
GE1D-PLOFPL-DMUO
GE1D-PLCFPL-DMUO
GE1D-PLPFPL-DMUO
GE1D-PLMFPL-DMUO
GE1D-PLF4PL-DMUO

Built to order. Call for lead time.

Add \$7 to substitute a PL (plug) for PT (90° Elbow Down DISS with check).

Substitute U with I for ISO color codes at no extra charge.

‡ Add \$30 and change 'GFS' to 'GFPS' to substitute Fairfield CAM adapter with push button style.



Double Gas Blocks - Medical Air

Replace (X) with selected adapter and (Y) with selected outlet. Shipped fully assembled with selected adapter style

Cost effective alternative to adding new medical gas outlets or moving outlets, suited for all gases.



G(X)D-PL(Y)PL-DMUA



By Adapter Style

Amico Horizontal Rail - A

Part#

GAD-PLDCPL-DMUA
 GAD-PLOFPL-DMUA
 GAD-PLCFPL-DMUA
 GAD-PLPFPL-DMUA
 GAD-PLMFPL-DMUA
 GAD-PLF4PL-DMUA

Amico Vertical Rail - V

Part#

GVD-PLDCPL-DMUA
 GVD-PLOFPL-DMUA
 GVD-PLCFPL-DMUA
 GVD-PLPFPL-DMUA
 GVD-PLMFPL-DMUA
 GVD-PLF4PL-DMUA

Direct To Wall - W

Part#

GWD-PLDCPL-DMUA
 GWD-PLOFPL-DMUA
 GWD-PLCFPL-DMUA
 GWD-PLPFPL-DMUA
 GWD-PLMFPL-DMUA
 GWD-PLF4PL-DMUA

Fairfield Rail - F

Part#

‡ GFD-PLDCPL-DMUA
 ‡ GFD-PLOFPL-DMUA
 ‡ GFD-PLCFPL-DMUA
 ‡ GFD-PLPFPL-DMUA
 ‡ GFD-PLMFPL-DMUA
 ‡ GFD-PLF4PL-DMUA

DIN 25mm Rail - E1

Part#

GE1D-PLDCPL-DMUA
 GE1D-PLOFPL-DMUA
 GE1D-PLCFPL-DMUA
 GE1D-PLPFPL-DMUA
 GE1D-PLMFPL-DMUA
 GE1D-PLF4PL-DMUA

Built to order. Call for lead time.

Add \$7 to substitute a PL (plug) for PT (90° Elbow Down DISS with check).

Substitute U with I for ISO color codes at no extra charge.

‡ Add \$30 and change 'GFS' to 'GFPS' to substitute Fairfield CAM adapter with push button style.

Double Gas Blocks - Vacuum

Replace (X) with selected adapter and (Y) with selected outlet. Shipped fully assembled with selected adapter style

Cost effective alternative to adding new vacuum outlets or moving outlets, suited for all gases.



G(X)D-PL(Y)PL-DMUV

Adapter Style (X)			Outlet Type (Y)		
A Amico Horizontal	F Fairfield	V Vertical	DC DISS	OF Ohmeda Compatible	CF Chemetron Compatible
E1 DIN 25mm	W Wall / Slide		PF Puritan Bennett Compatible	MF Medstar Compatible	F4 1/4" FNPT Compatible

By Adapter Style

Amico Horizontal Rail - A

Part#

GAD-PLDCPL-DMUV
GAD-PLOFPL-DMUV
GAD-PLCFPL-DMUV
GAD-PLPFPL-DMUV
GAD-PLMFPL-DMUV
GAD-PLF4PL-DMUV

Amico Vertical Rail - V

Part#

GVD-PLDCPL-DMUV
GVD-PLOFPL-DMUV
GVD-PLCFPL-DMUV
GVD-PLPFPL-DMUV
GVD-PLMFPL-DMUV
GVD-PLF4PL-DMUV

Direct To Wall - W

Part#

GWD-PLDCPL-DMUV
GWD-PLOFPL-DMUV
GWD-PLCFPL-DMUV
GWD-PLPFPL-DMUV
GWD-PLMFPL-DMUV
GWD-PLF4PL-DMUV

Fairfield Rail - F

Part#

‡ GFD-PLDCPL-DMUV
‡ GFD-PLOFPL-DMUV
‡ GFD-PLCFPL-DMUV
‡ GFD-PLPFPL-DMUV
‡ GFD-PLMFPL-DMUV
‡ GFD-PLF4PL-DMUV

DIN 25mm Rail - E1

Part#

GE1D-PLDCPL-DMUV
GE1D-PLOFPL-DMUV
GE1D-PLCFPL-DMUV
GE1D-PLPFPL-DMUV
GE1D-PLMFPL-DMUV
GE1D-PLF4PL-DMUV

Built to order. Call for lead time.

Add \$7 to substitute a PL (plug) for PT (90° Elbow Down DISS with check).

Substitute U with I for ISO color codes at no extra charge.

‡ Add \$30 and change 'GFS' to 'GFPS' to substitute Fairfield CAM adapter with push button style.


Triple Gas Blocks

Replace (X) with selected adapter and (Y) with selected outlet. Shipped fully assembled with selected adapter style

Cost effective alternative to adding new medical gas outlets or moving outlets, suited for all gases.



G(X)T-PL(Y)PL-DMU(Z)



Oxygen (By Adapter Style)

Fairfield Rail - F

Part#

- ‡ GFT-PLDCPL-DMUO
- ‡ GFT-PLOFPL-DMUO
- ‡ GFT-PLPFPL-DMUO
- ‡ GFT-PLF4PL-DMUO

Amico Vertical Rail - V

Part#

- GVT-PLDCPL-DMUO
- GVT-PLOFPL-DMUO
- GVT-PLPFPL-DMUO
- GVT-PLF4PL-DMUO

DIN 25mm Rail - E1

Part#

- GE1T-PLDCPL-DMUO
- GE1T-PLOFPL-DMUO
- GE1T-PLPFPL-DMUO
- GE1T-PLF4PL-DMUO

Direct To Wall - W

Part#

- GWT-PLDCPL-DMUO
- GWT-PLOFPL-DMUO
- GWT-PLPFPL-DMUO
- GWT-PLF4PL-DMUO

Medical Air (By Adapter Style)

Fairfield Rail - F

Part#

- ‡ GFT-PLDCPL-DMUA
- ‡ GFT-PLOFPL-DMUA
- ‡ GFT-PLPFPL-DMUA
- ‡ GFT-PLF4PL-DMUA

Amico Vertical Rail - V

Part#

- GVT-PLDCPL-DMUA
- GVT-PLOFPL-DMUA
- GVT-PLPFPL-DMUA
- GVT-PLF4PL-DMUA

DIN 25mm Rail - E1

Part#

- GE1T-PLDCPL-DMUA
- GE1T-PLOFPL-DMUA
- GE1T-PLPFPL-DMUA
- GE1T-PLF4PL-DMUA

Direct To Wall - W

Part#

- GWT-PLDCPL-DMUA
- GWT-PLOFPL-DMUA
- GWT-PLPFPL-DMUA
- GWT-PLF4PL-DMUA

Built to order. Call for lead time.

Add \$7 to substitute a PL (plug) for PT (90° Elbow Down DISS with check).

Substitute U with I for ISO color codes at no extra charge.

- ‡ Add \$30 and change 'GFS' to 'GFPS' to substitute Fairfield CAM adapter with push button style.


T Connector

Replace (X) with selected inlet, (Y) with selected outlet and (Z) with gas type.

Allows for the existing gas outlet to remain in use while adding a second gas channel.



GGB-(X)(Y)-DDCU(Z)



By Gas Type

Oxygen

Part#

GGB-DHDC-DDCUO
GGB-OMOF-DDCUO
GGB-CMCF-DDCUO
GGB-PMPF-DDCUO

Medical Air

Part#

GGB-DHDC-DDCUA
GGB-OMOF-DDCUA
GGB-CMCF-DDCUA
GGB-PMPF-DDCUA

Vacuum

Part#

GGB-DHDC-DDCUV
GGB-OMOF-DDCUV
GGB-CMCF-DDCUV
GGB-PMPF-DDCUV

Substitute U with I for ISO color codes at no extra charge.


Hoses available upon request. Please contact us for details.





Y Connector

Replace (X) with selected outlet and (Y) with gas type.

GGB-(X)-DCDCU(Y)

Connects to the medical gas source for dual gas distribution.



Outlet Type (X)				Gas Type (Y)		
				O	A	V
DISS Hand Tight by two DISS	Ohmeda Male by two DISS	Chemetron Male by two DISS	Puritan Bennett Male by two DISS	Oxygen	Medical Air	Vacuum
DH	OM	CM	PM			

By Gas Type

Oxygen

Part#

GGB-DH-DCDCUO
GGB-OM-DCDCUO
GGB-CM-DCDCUO
GGB-PM-DCDCUO

Medical Air

Part#

GGB-DH-DCDCUA
GGB-OM-DCDCUA
GGB-CM-DCDCUA
GGB-PM-DCDCUA

Vacuum

Part#

GGB-DH-DCDCUV
GGB-OM-DCDCUV
GGB-CM-DCDCUV
GGB-PM-DCDCUV


Substitute U with I for ISO color codes at no extra charge. Hoses available upon request. Please contact us for details.






Movable Horizon Gas Block

Replace (X) with selected outlet and (Y) with gas type.

GHZ-(X)06-U(Y)

Provides a moveable outlet on the Hill-Rom Horizon Headwall, includes 6' (1.83m) hose assembly.



Outlet Type (X)					Gas Type (Y)		
					O	A	V
DC DISS	Ohmeda Compatible	Chemetron Compatible	Puritan Bennett Compatible	Medstar Compatible	Oxygen	Medical Air	Vacuum

By Gas Type

Oxygen

Part#

GHZ-DC06-UO
GHZ-OF06-UO
GHZ-CF06-UO
GHZ-PF06-UO
GHZ-MF06-UO

Medical Air

Part#

GHZ-DC06-UA
GHZ-OF06-UA
GHZ-CF06-UA
GHZ-PF06-UA
GHZ-MF06-UA

Vacuum

Part#

GHZ-DC06-UV
GHZ-OF06-UV
GHZ-CF06-UV
GHZ-PF06-UV
GHZ-MF06-UV

Built to order. Call for lead time.

Substitute U with I for ISO color codes at no extra charge. Hoses at custom length available. Please contact us for details.

Amico Horizontal Adapter

ATA-8-32

Amico horizontal track adapter compatible with Amico and Modular horizontal rail system. Features (4) 8-32 threaded holes 1/2" apart for use on custom accessories applications.



Part#

ATA-8-32

Consult with manufacturer for maximum weight limit.

Amico Horizontal Socket Adapter

ATA-05

Amico horizontal track adapter compatible with Amico and Modular horizontal rail system. Features 1/2" socket for use with custom accessories.



Part#

ATA-05

Consult with manufacturer for maximum weight limit.

Amico Vertical Adapter

VTA-8-32

Amico vertical track adapter compatible with Amico and Modular vertical rail system. Features (2) 8-32 threaded holes 1/2" apart for use on custom accessories applications.



Part#

VTA-8-32

Consult with manufacturer for maximum weight limit.

Amico Vertical Socket Adapter

VTA-05

Amico vertical track adapter compatible with Amico and Modular vertical rail system. Features 1/2" socket for use with custom accessories.



Part#

VTA-05

Consult with manufacturer for maximum weight limit.

Fairfield Adapter

FTA-8-32

Amico Fairfield compatible adapter. Features (2) 8-32 threaded holes 1/2" apart for use on custom accessories applications.



Part#

FTA-8-32

Consult with manufacturer for maximum weight limit.

Fairfield Socket Adapter

FTA-05

Amico Fairfield compatible adapter with set screw lock down. Features 1/2" socket for use with custom accessories.



Part#

FTA-05

Consult with manufacturer for maximum weight limit.

Euro Track (DIN) 25mm Adapter

ETA-25-8-32

Amico DIN track adapter used with 10mm x 25mm European (DIN) rails with push button release mechanism. Features (2) 8-32 threaded holes 1/2" apart for use on custom accessories applications.



Part#

ETA-25-8-32

Consult with manufacturer for maximum weight limit.

Euro Track (DIN) 25mm Socket Adapter

ETA-25-05

Amico DIN track adapter used with 10mm x 25mm European (DIN) rails with push button release mechanism. Features 1/2" socket for use with custom accessories.



Part#

ETA-25-05

Consult with manufacturer for maximum weight limit.

MRS Vertical Monitor Channel Adapter

MHT-8-32

Amico monitor channel adapter which is compatible with horizontally or vertically mounted monitor channel. Features (2) 8-32 threaded holes 1/2" apart for use on custom accessories applications.



Part#

MHT-8-32

Consult with manufacturer for maximum weight limit.

MRS Vertical Monitor Channel Socket Adapter

MHT-05

Amico monitor channel adapter which is compatible with horizontally or vertically mounted monitor channel. Features 1/2" socket for use with custom accessories.



Part#

MHT-05

Consult with manufacturer for maximum weight limit.

Hill-Rom Vertical Adapter

GTA-8-32

Amico track adapter compatible with Hill-Rom vertical rail system. Features (2) 8-32 threaded holes 1/2" apart for use on custom accessories applications.



Part#

GTA-8-32

Consult with manufacturer for maximum weight limit.

Hill-Rom Vertical Socket Adapter

GTA-05

Amico track adapter compatible with Hill-Rom vertical rail system. Features 1/2" socket for use with custom accessories.



Part#

GTA-05

Consult with manufacturer for maximum weight limit.

HSI Vertical Rail Adapter

STA-8-32

Amico vertical track adapter compatible with HSI vertical rail system. Features (2) 8-32 threaded holes 1/2" apart for use on custom accessories applications.



Part#

* STA-8-32

* Long lead time product.

Consult with manufacturer for maximum weight limit.

HSI Vertical Rail Socket Adapter

STA-05

Amico vertical track adapter compatible with HSI vertical rail system. Features 1/2" socket for use with custom accessories.



Part#

* STA-05

* Long lead time product.

Consult with manufacturer for maximum weight limit.

Ohmeda Dovetail Adapter

OHM-8-32

Amico dovetail adapter compatible with Datex-Ohmeda/GE dovetail rail systems. Features (2) 8-32 threaded holes 1/2" apart for use on custom accessories applications.



Part#

OHM-8-32

Consult with manufacturer for maximum weight limit.

Ohmeda Dovetail Socket Adapter

OHM-05

Amico dovetail adapter compatible with Datex-Ohmeda/GE dovetail rail systems. Features 1/2" socket for use with custom accessories.



Part#

OHM-05

Consult with manufacturer for maximum weight limit.

Direct Wall Adapter

DWA-8-32

Amico direct wall adapter used for mounting directly to wall, or into an Amico vacuum bottle slide. Features (2) 8-32 threaded holes 1/2" apart for use on custom accessories applications and (2) 0.156 TYP holes 3.28" apart for direct wall mount.



Part#

* DWA-8-32

* Long lead time product.

Consult with manufacturer for maximum weight limit.

Direct Wall Socket Adapter

DWA-05

Amico direct wall adapter used for mounting directly to wall, or into an Amico vacuum bottle slide. Features 1/2" socket for use with custom accessories and (2) 0.156 TYP holes 3.28" apart for direct wall mount.



Part#

* DWA-05

* Long lead time product.

Consult with manufacturer for maximum weight limit.

Adjustable Pole Clamp

PC1-8-32

Amico adjustable C-clamp, compatible with poles at 0.5" (1.27cm) diameter to 1.5" (3.8cm). Features (4) 8-32 threaded holes 1/2" apart for use on custom accessories applications.



Part#

* PC1-8-32

* Long lead time product.

Consult with manufacturer for maximum weight limit.

C Clamp for 1.25" Post

PC-125-8-32

Amico 1.25" C-clamp, compatible with poles at 1.25"(3.18cm) diameter. Features (4) 8-32 threaded holes 1/2" apart for use on custom accessories applications.



Part#

PC-125-8-32

Consult with manufacturer for maximum weight limit.

C Clamp for 1.5" Post

PC-150-8-32

Amico 1.5" C-clamp, compatible with poles at 1.5"(3.81cm) diameter. Features (4) 8-32 threaded holes 1/2" apart for use on custom accessories applications.



Part#

* PC-150-8-32

* Long lead time product.

Consult with manufacturer for maximum weight limit.

C Clamp for 2" Post

PC-2-8-32

Amico 2" C-clamp, compatible with poles at 2"(5.08cm) diameter. Features (4) 8-32 threaded holes 1/2" apart for use on custom accessories applications.



Part#

* PC-2-8-32

* Long lead time product.

Consult with manufacturer for maximum weight limit.

Amico Horizontal Rail System

Replace (X) with length by inch
Shipped fully assembled
with end caps

ARS-(XXX)

Designed with negative profile allows for ease of maintenance and infection control available in a wide variety of lengths. Sold by foot length.



Part#

ARS-012
ARS-024
ARS-036
ARS-048
ARS-060

Part#

ARS-072
ARS-084
ARS-096
ARS-108
ARS-120

Max. load: 15 lb. ft. with typical mounting arrangement (2x 1/4" TEK fastener through the ARS/LRS/FRS bracket into gauge 16 steel backing every 16". The weight limit applies between every 3 studs used (i.e. 3 stud gap before another load can be applied on the rail).

Fasteners required to mount to wall are not provided.

Extra shipping charges apply to rail system over 72" as special carrier is required.

Light Duty Fairfield Rail System

Replace (X) with length by inch
Shipped fully assembled
with end caps

LRS-(XXX)

Simple, functional and economical solution designed with a slim profile ideal for adding to existing equipment to match variable needs within the patient care environment. Sold by foot length.



Part#

LRS-012
LRS-024
LRS-036
LRS-048
LRS-060

Part#

LRS-072
LRS-084
LRS-096
LRS-108
LRS-120

Max. load: 15 lb. ft. with typical mounting arrangement (2x 1/4" TEK fastener through the ARS/LRS/FRS bracket into gauge 16 steel backing every 16". The weight limit applies between every 3 studs used (i.e. 3 stud gap before another load can be applied on the rail).

Fasteners required to mount to wall are not provided.

Extra shipping charges apply to rail system over 72" as special carrier is required.

Fairfield Rail System

Replace (X) with length by inch
Shipped fully assembled
with end caps

FRS-(XXX)

Extended profile for support on horizontal applications. Available to be fitted directly to the wall surface or paired with mounting hardware.



Part#

FRS-012
FRS-024
FRS-036
FRS-048
FRS-060

Part#

FRS-072
FRS-084
FRS-096
FRS-108
FRS-120

Max. load: 15 lb. ft. with typical mounting arrangement (2x 1/4" TEK fastener through the ARS/LRS/FRS bracket into gauge 16 steel backing every 16". The weight limit applies between every 3 studs used (i.e. 3 stud gap before another load can be applied on the rail).

Fasteners required to mount to wall are not provided.

Extra shipping charges apply to rail system over 72" as special carrier is required.

Extended Wall Mounting Bracket for Light Duty/Fairfield Rail System

FRS-E-WALL-BRKT

Extended wall mounting bracket for Light Duty/Fairfield Rail System, 2 brackets minimum per rail. Recommended a minimum of 1 bracket for every 16" length.



Part#

FRS-E-WALL-BRKT

Mounting hardware not included.

Direct Wall Mounting Bracket for Light Duty/Fairfield Rail System

FRS-D-WALL-BRKT

Direct wall mounting bracket for Light Duty/Fairfield rail system, 2 brackets minimum per rail. Recommended a minimum of 1 bracket for every 16" length.



Part#

FRS-D-WALL-BRKT

Mounting hardware not included.

Light Duty Rail System (Conversion)

Replace (X) with length by inch
Shipped fully assembled with
end caps and adapters

LRS-(XXX)-A

A - Amico Horizontal Rail

Simple and economical solution to convert Amico horizontal rail systems into light duty Fairfield rail system. Sold by foot length.



Part#

LRS-012-A
LRS-024-A
LRS-036-A
LRS-048-A
LRS-060-A

Part#

LRS-072-A
LRS-084-A
LRS-096-A
LRS-108-A
LRS-120-A

Extra shipping charges apply to rail system over 72" as special carrier is required.

Amico Vertical Rail System

Replace (X) with length by foot
Shipped fully assembled
with end caps

VRS-(XXX)

Amico Vertical Rail System.



Part#

VRS-012
VRS-024
VRS-036
VRS-048
VRS-060

Part#

VRS-072
VRS-084
VRS-096
VRS-108
VRS-120

Fasteners required to mount to wall are not provided.

Extra shipping charges apply to rail system over 72" as special carrier is required.

Amico Horizontal Rail System on Adapter for Vertical Flatwall

Replace (X) with length
Shipped fully assembled with
end caps and adapters

ARS-(XXX)-V

Perfect for adding ARS horizontal rail versatility to existing vertical systems, allows users to maximize unused space.



Part#	Mounting Application
ARS-009-V	9" Amico Clinical Boom Console Spanning 2 Vertical Rails (7.3" c-c)
ARS-012-V	12" Amico Power Column/Chase Spanning 2 Vertical Rails (10.3" c-c)
ARS-012-V	12.5" Amico Clinical Boom Console Spanning 2 Vertical Rails (10.8" c-c)
ARS-014-V	14" Amico Clinical Boom Console Spanning 2 Vertical Rails (12.3" c-c)
ARS-015-V	14" Amico Recessed Flatwall Spanning 2 Vertical Rails (13.9" c-c)
ARS-018-V	18" Amico Clinical Boom Console Spanning 2 Vertical Rails (16.3" c-c)
ARS-019-V	19" Amico Surface Flatwall Spanning 2 Vertical Rails (17.8" c-c)
ARS-020-V	19" Amico Recessed Flatwall Spanning 2 Vertical Rails (18.9" c-c)
ARS-024-V	24" Amico Surface Flatwall Spanning 2 Vertical Rails (22.3" c-c)
ARS-028-V	28" Amico Surface Flatwall Spanning 2 Vertical Rails (26.8" c-c)

Max. load: 10 lb. ft. with typical mounting arrangement (2x 1/4" TEK fastener through the ARS/LRS/FRS bracket into gauge 16 steel backing every 16". The weight limit applies between every 3 studs used (i.e. 3 stud gap before another load can be applied on the rail).)

Light Duty Fairfield Rail System on Adapter for Vertical Flatwall

Replace (X) with length
Shipped fully assembled with
end caps and adapters

LRS-(XXX)-V

Perfect for adding LRS horizontal rail versatility to existing vertical systems, allows users to maximize unused space.



Part#	Mounting Application
LRS-011-V	9" Amico Clinical Boom Console Spanning 2 Vertical Rails (7.3" c-c)
LRS-012-V	12" Amico Power Column/Chase Spanning 2 Vertical Rails (10.3" c-c)
LRS-014-V	12.5" Amico Clinical Boom Console Spanning 2 Vertical Rails (10.8" c-c)
LRS-016-V	14" Amico Clinical Boom Console Spanning 2 Vertical Rails (12.3" c-c)
LRS-016-V	14" Amico Recessed Flatwall Spanning 2 Vertical Rails (13.9" c-c)
LRS-020-V	18" Amico Clinical Boom Console Spanning 2 Vertical Rails (16.3" c-c)
LRS-020-V	19" Amico Surface Flatwall Spanning 2 Vertical Rails (17.8" c-c)
LRS-021-V	19" Amico Recessed Flatwall Spanning 2 Vertical Rails (18.9" c-c)
LRS-024-V	24" Amico Surface Flatwall Spanning 2 Vertical Rails (22.3" c-c)
LRS-029-V	28" Amico Surface Flatwall Spanning 2 Vertical Rails (26.8" c-c)

Max. load: 10 lb. ft. with typical mounting arrangement (2x 1/4" TEK fastener through the ARS/LRS/FRS bracket into gauge 16 steel backing every 16". The weight limit applies between every 3 studs used (i.e. 3 stud gap before another load can be applied on the rail).)

Fairfield Rail System on Adapter for Vertical Flatwall

Replace (X) with length
Shipped fully assembled with
end caps and adapters
FRS-(XXX)-V

Cost effective solution to add horizontal Fairfield rail system on Amico vertical systems.



Part#	Mounting Application
FRS-011-V	9" Amico Clinical Boom Console Spanning 2 Vertical Rails (7.3" c-c)
FRS-012-V	12" Amico Power Column/Chase Spanning 2 Vertical Rails (10.3" c-c)
FRS-014-V	12.5" Amico Clinical Boom Console Spanning 2 Vertical Rails (10.8" c-c)
FRS-016-V	14" Amico Clinical Boom Console Spanning 2 Vertical Rails (12.3" c-c)
FRS-016-V	14" Amico Recessed Flatwall Spanning 2 Vertical Rails (13.9" c-c)
FRS-020-V	18" Amico Clinical Boom Console Spanning 2 Vertical Rails (16.3" c-c)
FRS-020-V	19" Amico Surface Flatwall Spanning 2 Vertical Rails (17.8" c-c)
FRS-021-V	19" Amico Recessed Flatwall Spanning 2 Vertical Rails (18.9" c-c)
FRS-024-V	24" Amico Surface Flatwall Spanning 2 Vertical Rails (22.3" c-c)
FRS-029-V	28" Amico Surface Flatwall Spanning 2 Vertical Rails (26.8" c-c)

Max. load: 10 lb. ft. with typical mounting arrangement (2x 1/4" TEK fastener through the ARS/LRS/FRS bracket into gauge 16 steel backing every 16". The weight limit applies between every 3 studs used (i.e. 3 stud gap before another load can be applied on the rail.)

Monitor Channel System

Replace (X) with length
Shipped fully assembled
MRS-(XXX)

Amico Monitor Rail System compatible with GCX monitor channel. Flanged.



Part#

MRS-024
MRS-034

Part#

MRS-XXX

80 lb. ft. (36.3 kgs) max. load at rail.

Custom Note: Custom length available for up to 120" and discount is unavailable for custom product. Sold by the next available longer foot length - LIST PRICE ONLY.

Monitor Channel System - Flangeless White 10"

MRS-010

Amico Monitor Rail System, powder coated white and flangeless, ideal for mounting small items like CPUs and Flush LCD monitors.



Part#

MRS-010

40 lb. ft. (18.2 kgs) max. load at rail.

Monitor Channel System - Flangeless

Replace (X) with length
Shipped fully assembled
MRS-FL-(XXX)

Amico Monitor Rail System flangeless and compatible with GCX monitor channel.



Part#

MRS-FL-019
MRS-FL-024

Part#

MRS-FL-034

40 lb. ft. (18.2 kgs) max. load at rail.

Monitor Channel for Draeger

DR-MRS-FL-006

6" accessory track for Draeger Apollo, Primus, Primus IE, Julian and Narkomed anesthesia.



Part#

* DR-MRS-FL-006

* Long lead time product.

Monitor Channel System - Seismic

Replace (X) with length
Shipped fully assembled

Amico Monitor Rail System compatible with GCX monitor channel - Seismic applications.



MRS-S(XXX)

Part#

MRS-S019
MRS-S024
MRS-S034

Part#

MRS-S048
MRS-SXXX

80 lb. ft. (36.3 kg) max. load for MRS-S019, 60 lb. ft. (27.2 kg) max. load for MRS-S034.
Custom Note: Custom length available for up to 48" and discount is unavailable for custom product. Sold by the next available longer length - LIST PRICE ONLY.

Monitor Channel Cover

Replace (X) with length

MRS-COVER-(XXX)

Ideal for concealing cables on any Amico monitor rail system (MRS) track. Flangeless cover, in white.



Part#

MRS-COVER-010
MRS-COVER-019
MRS-COVER-024

Part#


MRS-COVER-034
MRS-COVER-048

Monitor Channel Cover - Flanged

Replace (XX) with length, (YYY) with Amico arm type, (Z) with colour

MRS-CVR-(XX)-(YYY)-(Z)

This flanged monitor channel cover offers ample space for cable management to create an organized and safe work environment. Standard in white and black.



BLACK - for AHM arm

Part#

MRS-CVR-19-AHM-B
MRS-CVR-34-AHM-B
MRS-CVR-48-AHM-B

BLACK - for LMA/SSM arm

Part#

MRS-CVR-19-LMA-B
MRS-CVR-34-LMA-B
MRS-CVR-48-LMA-B

WHITE - for AHM arm

Part#

MRS-CVR-19-AHM-W
MRS-CVR-34-AHM-W
MRS-CVR-48-AHM-W

WHITE - for LMA/SSM arm

Part#


MRS-CVR-19-LMA-W
MRS-CVR-34-LMA-W
MRS-CVR-48-LMA-W

Monitor Channel for Horizontal Headwall

Replace (X) with selected adapter
Shipped fully assembled with
selected adapter style

MRS-020-(X)

Easily allows facilities to alleviate the need for wall backing and costly installation when mounting to a horizontal headwall system.



F - Fairfield Rail

A - (18" c-c) Spans 2 Amico Horizontal Rails 18" Apart

A75 - Amico (7.5" c-c) Spans 2 Horizontal Rails

ARA - (15-18" c-c) Spans 2 Amico Horizontal Rails with Standoff

Part#	Mounting Application
* MRS-020-F	Spans two (2) Fairfield rails 15"-18" apart
MRS-020-A	Spans two (2) Amico horizontal rails 15"-18" apart, or span across the 1st to 3rd rail of a Amico Majestic 2-tier headwall
MRS-020-A75	Spans two (2) Amico horizontal rails 7.5" apart, or span across the 1st to the 2nd rail of a Amico Majestic 1-tier Headwall
* MRS-020-ARA	Spans from a surface mount Amico horizontal rail to a Amico headwall horizontal rail, 15" to 18" apart

60 lb. ft. (27.2 kgs) max. load. Please always check the center-to-center distance between the two mounting points.

* Long lead time product.

Stud Spanning Kit for MRS

MRS-034-STSP

Stud Spanning Kit for 34" MRS rail (FRS mount)



Part#

MRS-034-STSP

Monitor Channel with Bumper

MRS-034-SMS

Monitor Channel with bumper for headwall mounting for Majestic Series Headwalls.



Part#

* MRS-034-SMS

* Long lead time product.

Monitor Channel on Hill-Rom Headwall

HR-MRS-006-MNT

6" Monitor Channel for headwall mounting to horizontal rail of Hill-Rom Headwall.



Part#

* HR-MRS-006-MNT

* Long lead time product.

Monitor Channel Stopper

CHAN-STOPPER

Monitor Channel Stopper



Part# _____

* CHAN-STOPPER

* Long lead time product.

Mobile Amico Horizontal Rail

Replace (X) with length by inch
Shipped fully assembled with
end caps and adapters

ARS-RS042-(XXX)

Mobile Amico Horizontal Rail mounted back-to-back on a rolling stand, offering flexibility in gas outlet when pair up with gas blocks.



Part# _____

* ARS-RS042-012

* ARS-RS042-024

* Long lead time product.

Mobile Fairfield Rail

Replace (X) with length by inch
Shipped fully assembled with
end caps and adapters

FRS-RS042-(XXX)

Mobile Fairfield Rail mounted back-to-back on a rolling stand, offering flexibility in gas outlet when pair up with gas blocks.



Part# _____

* FRS-RS042-012


* FRS-RS042-024

* Long lead time product.

Monitor Roll Stand

Replace (X) with selected counterweight and (Y) with selected head. Shipped fully assembled with the selected counterweight and head style.

42"(106.7cm) patient monitor roll stand with 5-wheel base.



By Counterweight Style

42" Tall Roll Stand with white post

RS5-42(Y)

Part#

RS5-420
RS5-422
RS5-423
RS5-424
RS5-427
RS5-428
RS5-429
RS5-42W

42" Tall Roll Stand with white post & counterweight

RS6-42(Y)

Part#

RS6-420
RS6-422
RS6-423
RS6-424
RS6-427
RS6-428
RS6-429
RS6-42W

42" Tall Roll Stand with white post, base & leg counterweight

RS7-42(Y)

Part#

RS7-420
RS7-422
RS7-423
RS7-424
RS7-427
RS7-428
RS7-429
RS7-42W


Monitor Roll Stand Package

Replace (X) with selected counterweight and (Y) with selected head. Shipped fully assembled with the selected counterweight and head style.


42"(106.7cm) patient monitor roll stand with 5-wheel base, includes basket and handle.

RS(X)-42(Y)-PKG

Counterweight Style (X)



Head Type (Y)



By Counterweight Style

42" Tall Roll Stand package with white post

RS5-42(Y)-PKG

Part#

RS5-420-PKG
RS5-422-PKG
RS5-423-PKG
RS5-424-PKG
RS5-427-PKG
RS5-428-PKG
RS5-429-PKG
RS5-42W-PKG

42" Tall Roll Stand package with white post & counterweight

RS6-42(Y)-PKG

Part#

RS6-420-PKG
RS6-422-PKG
RS6-423-PKG
RS6-424-PKG
RS6-427-PKG
RS6-428-PKG
RS6-429-PKG
RS6-42W-PKG

42" Tall Roll Stand package with white post, base & leg counterweight

RS7-42(Y)-PKG

Part#


RS7-420-PKG
RS7-422-PKG
RS7-423-PKG
RS7-424-PKG
RS7-427-PKG
RS7-428-PKG
RS7-429-PKG
RS7-42W-PKG

LCD Roll Stand

Replace (X) with selected counterweight. Shipped fully assembled with selected counterweight style

RS(X)-LCD

42"(106.7cm) LCD/flat panel (75/100mm VESA) roll stand with 5-wheel base.



Part#

RS5-LCD
RS6-LCD

Part#


RS7-LCD

LCD Roll Stand Package

Replace (X) with selected counterweight. Shipped fully assembled with selected counterweight style

RS(X)-LCD-PKG

42"(106.7cm) LCD/flat panel (75/100mm VESA) roll stand with 5-wheel base, includes basket and handle.



Part#

RS5-LCD-PKG
RS6-LCD-PKG

Part#

RS7-LCD-PKG